



# Agilent 6140 UHPLC-MS

-Scott Virgil, California Institute of Technology May, 2016




The use of UHPLC technology in combination with fast Electrospray Mass Spectrometry has advanced the speed and efficiency of research in demanding fields of Organic Synthesis and Catalysis.

In the above diagram, the key features of the Agilent 6140 UHPLC-MS instrument are shown:

1. The **1290 Binary Pump** module delivers water and acetonitrile at pressures up to 1200 bar. It also includes a degasser, an auto-purge and a “Jet Weaver” microfluidic.
2. The **1290 Well-Plate Autosampler** accepts standard HPLC vial trays and 96-well plates.
3. The sample flows through a 2  $\mu\text{m}$  filter to a 2.1 x 50 mm Eclipse Plus 1.8  $\mu\text{m}$  particle size C18 column (part No. 959757-902) at a flow rate of 1 mL/minute and 700-800 bar.
4. The **1290 Diode Array Detector (DAD)** with a 60 mm path length microfluidic flowcell detects the absorbance a 6x improvement in sensitivity over standard 10 mm flowcells.
5. The **6140 MS Single Quad MS detector** with standard ESI source is specifically designed to provide fast data acquisition on a mass range of 50-1350 m/e range.


## Comparison of 1 ng injection on 1100 LCMS vs. 1290 UHPLC-MS:

- 2-4x Improvement in Resolution, 10x Improvement in Sensitivity



# UHPLC-MS System Operation with Agilent Chemstation and Easy-Access

The Agilent UHPLC-MS is operated using an **Online Chemstation** window that runs the instrument and an **Offline Chemstation** window that is used for *data analysis*. In addition, **Agilent Easy-Access** software is used to efficiently handle multiple users in a walk-up format.



## I. How to submit samples to the Easy Access queue.

- Hit “Submit New Sample” and enter your login account and password.
- Select the number of samples to submit and enter the **Sample ID(s)**.

MassHunter Easy-Access - Sample Data Input - SCV

User Login		Sample Data Input			Sample Loading
SAMPLE ID	Method	Inj. Vol.	UV (nm)		
1 SCV-III-084	PES10BX	1.00	254		
2 SCV-III-085	PES4AX_20-95	1.20	272		
3 SCV-III-087 "crude"	NES7CX	1.00	254		
4 SCV-III-084	NES4AX_MSOpen0min	1.00	254		

Examples - **Correct**  
**X - Don't use . or \* in Sample Name**  
**X - Don't duplicate Sample Name**

☒ Automatically copy down column

Back Load Samples Cancel

• Choose from the set of available Methods. Either **Positive** (PES) or **Negative** (NES) detection can be used in LCMS electrospray. However, positive electrospray ionizes a wider range of compounds. The methods have run times of 4, 7 or 10 minute lengths. The “AX” methods have a **100-800 amu** mass range. For compounds that are above 190 amu, the “BX” has **190-800 amu** (which eliminates low MW impurities from the MS trace) and “CX” runs from **190-1350 amu**. If no other annotations are included, the gradient proceeds from 5% to 95% acetonitrile during the run time and the MS detector is bypassed until  $t = 0.5$  min in order to minimize salts entering the MSD. For polar compounds, the \_MSOpen0min can be used to open the MS at  $t = 0$  min. Both the default Inj. Vol. and the  $UV_A$  wavelength may also be changed.

PES10BX	Pos. 10 min run, 5-95% ACN, 190-800 amu
PES4AX_20-95	Pos. 4 min run, 20-95% ACN, 100-800 amu
NES7CX_MSOpen0min	Neg. 7 min run, 5-95% ACN, 190-1350 amu, MS valve open 0 min

## UHPLC-MS System Operation (continued)

### II. Common Questions on Samples:

- **How much sample do I need?** The LCMS easily detects analytes at the 1 ng level which corresponds to a 1  $\mu$ L injection from a 0.001 mg/mL sample. However, most samples are about 0.1 mg/mL and give UV peaks  $\sim$  500 mAu. Samples that are at or above 1 mg/mL are unnecessarily concentrated and will give poor data and diminish the life of the instrument.
- **What solvents are OK?** Acetonitrile is available for making up samples at the instrument. Most other standard organic solvents and water are suitable for sample submission. Chloroform and dichloromethane can affect the peak shapes if used at higher than 1  $\mu$ L injection volumes.
- **How much solution do I need?** The WP autosampler needs about 0.5 mL in a hplc vial as it draws from about 1/4 inch above the bottom. Vial inserts are available for smaller sample sizes.
- **Do I need to filter my sample?** No, unless the sample is obviously heterogeneous.
- **Can I submit “crude reaction mixtures”?** Yes. For reaction monitoring, usually a tip of a pipette diluted to 1 mL will provide a suitable sample.
- **What about biological samples?** Well, the worst thing for a C18 column to see is biological material including membrane proteins, phospholipids and biopolymers. These materials need to be removed from LCMS samples prior to submission or they will clog the column and degrade the peak shape of analytes.

#### Further notes on sample submission:

- When submitting more than 3 samples at a time, select the “Night Queue” to allow other samples to run first.
- For samples that are highly lipophilic, 20-95, 40-95 or 60-95 gradients can be used to improve separation and allow shorter run times.

### III. Data Analysis

1. **Auto-Print From Easy Access** – check the box “Print Report” when submitting samples.


2. **Snapshot** – while a sample is running, Select “Snapshot” in the Offline Chemstation.

3. **Load Signal** to access previously run Data – then select file on C:/ or D:/ (old data) drive. If only one signal chromatogram is desired, it may be selected individually when loading.

### IV. Tools for Data Analysis


When a data file is loaded, it will usually consist of three DAD chromatograms plus one MSD trace. They can be viewed individually or all together by selecting “All Loaded Signals”. In addition, five buttons above the data provide different tools for Integration, Signal selection, etc.

The “Spectrum” button allows one to view the Mass Spectrum or UV spectrum of a selected peak.




## UHPLC-MS Data Analysis (continued)

Shown below, the dragging button delivers a spectrum for the region that the user drags the mouse over. The UV spectrum is created by subtracting the red and blue UV absorbances from the dragged region in order to compensate for solvent absorbance (the acetic acid used in the aqueous mobile phase absorbs in the 210 – 230 nm range).




**Extract Ions** – used to search for a mass or to evaluate the MS purity of a peak:




## UHPLC-MS Data Analysis (continued)

**Modifying Data Analysis Preferences** – The first thing to know is how to get back to the default “LCMS\_ANALYSIS” Method since previous users often change settings.



### Methods for Integration:

1. Using the Integration button, a series of tools appear for manual integration.
2. Using “Integration Events”, the integration settings may be adjusted:



### Integration Events: showing how to raise the Area reject:

**Hit Green Check Mark when done.**

Method Manual Events: ☐

Initial Events For All Signals:

Integration Events	Value
Tangent Skim Mode	Standard
Tail Peak Skim Height Ratio	0.00
Front Peak Skim Height Ratio	0.00
Skim Valley Ratio	20.00
Baseline Correction	Advanced
Peak to Valley Ratio	500.00

Specific Events For Signal: MSD1 TIC Specific

Time	Integration Events	Value
Initial	Slope Sensitivity	100000
Initial	Peak Width	0.07
Initial	Area Reject	300000
Initial	Height Reject	50000
Initial	Shoulders	OFF
0.000	Integration	OFF
0.400	Integration	ON
10.000	Integration	OFF

**Select "Integrate" to re-integrate.**

MSD1 TIC, MS File (C:\CHEM32\1\DATA\SVIRGIL\06-16\0622TEST.D) ES-API, Pos, Scan, Frag: 70, 4

#	Time	Area	Height	Width	Area%	Symmetry
1	3.51	2540690	613036.8	0.0675	15.838	0.891
2	3.637	3524030	756434.4	0.0735	21.968	0.73
3	3.957	3072087	558283	0.0833	1.915	0.733
4	4.491	5182688	847486.1	0.0904	32.308	0.462
5	4.937	352658.2	66743.7	0.076	2.200	0.577
6	5.469	3785635.3	833323.9	0.0721	23.599	0.734
7	5.768	348297.2	85184.3	0.062	2.171	0.958



## UHPLC-MS Data Analysis (continued)

**Printing Report:** The “Specify Report” Window is used to select the format and destination of the printout. When re-loading the LCMSANALYSIS method, it is set to the printer. One can preview by selecting “Screen” or print to file formats such as pdf. When satisfied, Select **Report:Print Report**.

